GB IE

BODY ANALYSER SCALE

Operation and Safety Notes

FI

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SE

VÅG SOM MÄTER KROPPSFETTET

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




B

Introduction

Proper use.....	Page 6
Parts list.....	Page 6
Included items	Page 6
Technical information	Page 7

Safety advice	Page 7
----------------------------	---------------

Before getting started

Inserting / replacing the batteries.....	Page 8
Mounting the control unit.....	Page 8

Preparing for use

Setting the time.....	Page 9
Setting the alarm	Page 9
Setting the unit of weight	Page 9
Entering personal details	Page 9

How to use

Measuring weight	Page 10
Measuring weight, calorie requirement, body fat, body water, bone and muscle percentages.....	Page 10
Retrieving last recorded values	Page 11

What the interpreted readings mean.....	Page 11
--	----------------

Trouble-shooting	Page 11
-------------------------------	----------------

Cleaning	Page 12
-----------------------	----------------

Disposal.....	Page 12
----------------------	----------------

Body Analyser Scale

● Introduction

Familiarise yourself with the product before using it for the first time. In addition please refer carefully to the operating instructions below and the safety advice. Use the product only as described and for the indicated purpose. Keep these instructions in a safe place. If you pass the product on to anyone else, please ensure that you also pass on all the documentation.

● Proper use

These scales are intended for measuring the weight of people. The scales have a weighing range of up to 180 kg. The scales are also suitable for measuring body fat, body water, muscle and bone mass. Furthermore, it is possible to determine your daily calorie requirement by entering your gender, age, height and weight. Any other use or modification of the device constitute improper use and may result in injury and/or damage to the device. The manufacturer is not liable for damage caused by improper use. The device is not intended for commercial use.

● Parts list

Scales

- 1 Inlet for control unit
- 2 Electrodes
- 3 Platform
- 4 Metal contact
- 5 Infrared transmitter
- 6 Battery compartment

Removable control unit

- 7 Time/confirmation button
- 8 Alarm setting/count-up button
- 9 Button to set the weight units/count-down button
- 10 Start button
- 11 LC display
- 12 Mounting device

- 13 Metal contact
- 14 Infrared receiver
- 15 Battery compartment lid clip
- 16 Battery compartment

Pictograms in the LC display ¹¹

- 17 **AM** – Symbol for “mornings”
- 18 – Alarm activated/alarm settings
- 19 **PM** – Symbol for “afternoons”
- 20 **st:lb** – Weight in “stones:pounds”
- lb** – Weight in “pounds”
- 21 **kg** – Weight in kilograms
- 22 **KCAL** – Recommended daily calorie requirement
- 23 ‘ ” – Height in feet and inches
- 24 **cm** – Height in cm
- 25 **%** – Percentages (of body fat, muscles, body water, bone)
- 26 – Body water percentage
- 27 – Bone mass
- 28 – Body fat percentage
- 29 – Age group between 7-17 years of age (children mode)
- 30 – Muscle mass
- 31 – Age group between 18-99 years of age (adult mode)
- 32 – Does a lot of sport
- 33 – Physically demanding work
- 34 – Woman
- 35 – Standing job
- 36 – Man
- 37 – Sedentary work
- 38 – Does no sport

● Included items

Check that all the items are present and that the device is in perfect condition immediately after unpacking.

- 1 x Body fat scales
- 1 x Control unit
- 4 x Battery
- 2 x Screw
- 2 x Dowel
- 1 x Operating instructions

● Technical information

Battery:	2 x 1.5 V --- battery (AAA) Micro (transmitter) 2 x 1.5 V --- battery (AAA) Micro (receiver)
Platform:	31 x 31 cm, safety glass
Load capacity:	max. 180 kg (400 lbs)
Weight increment:	100 g (0.2 lb)
Minimum starting weight:	3 kg (6.6 lb)
Calorie increment:	1 kcal
Memory capacity:	12 + guest (total: 13)
Gender:	male / female
Age range:	7-99 years of age
Activity level:	5 (recumbent, sedentary, standing, physical effort, sports)
Height range:	100-250 cm (3'30 ½" - 8'02 ½")
Operating temperature:	10-40 °C
Measuring range:	
Body fat:	4.0-60.0 %
Body water:	27.5-66.0 %
Muscles:	20.0-56.0 %
Bones:	2.0-20.0 %

Safety advice

This safety advice is an important part of these operating instructions and must be observed every time electrical devices are used. The advice and instructions will help you to avoid risk of loss of life, accidents and damage. Therefore you must read carefully through the following information concerning safety and proper use.

- 
⚠ WARNING! RISK OF LOSS OF LIFE OR ACCIDENT TO INFANTS AND CHILDREN! Never leave children unsupervised with the packaging materials or the device. Otherwise there is a risk to life from suffocation. Always keep children away from the device.

- Children or persons who lack the knowledge or experience to use the device or whose physical, sensory or intellectual capacities are limited must never be allowed to use the device without supervision or instruction by a person responsible for their safety. Children must never be allowed to play with the device. The product is not a toy.

⚠ WARNING! DANGER TO LIFE! Never use the product if you have a pacemaker. Never use the product near people that have a pacemaker or use other medical devices.

⚠ CAUTION! DANGER OF INJURY! Do not load the scales beyond 180 kg. The safety glass could break and cause cut injuries

⚠ CAUTION! DANGER OF INJURY!

Do not step on to the scales with wet feet and do not use the scales if the platform is moist. You might slip.

- Protect your scales from impact, moisture, dust, chemicals and sudden changes in temperature and keep them away from sources of heat (ovens, heating radiators).
- Do not expose the device to:
 - extreme temperatures,
 - strong vibrations,
 - high mechanical loads,
 - direct solar radiation,
 - moisture.

Failure to observe this advice may result in damage to the device.

- Check the device for damage before bringing it into use. If the device is damaged do not use it.
- Remove the battery if the device breaks down or is defective. Otherwise the battery is more likely to leak.
- Never open the device housing. This could result in personal injury or damage to the device. Always have any repairs carried out at the service centre or by an electrical equipment repair specialist.
- The product is not suitable for pregnant women.

Warning! Dangers arising from the use of batteries

- ⚠ WARNING! DANGER TO LIFE!** Batteries are not intended to be in the hands of children. Do not leave batteries lying around. They could be swallowed by children or pets. In the event of swallowing a battery, contact a doctor immediately!

CAUTION! DANGER OF EXPLOSION! Never recharge single-use batteries. Do not short-circuit or open them. Overheating, fire or destruction of the batteries might result. Never throw batteries into fire or water. The batteries might explode.

- Avoid extreme environmental conditions and temperatures as they may affect batteries e.g. above radiators. Otherwise the batteries are more likely to leak.
- Exhausted or damaged batteries can cause caustic burns upon contact with the skin. To avoid this, always wear suitable protective gloves in these circumstances!
- Avoid contact with the skin, eyes and mucous membranes if fluid is leaking from a battery. If you come into contact with corrosive battery fluids, flush the affected area thoroughly with clean water and contact a doctor immediately.
- Clean any light dirt deposits on the contacts of the batteries or on the contacts of the device with a dry, fluff-free cloth before insertion.
- Remove the battery if you are not going to use the device for an extended period.
- When inserting the batteries, ensure the correct polarity! This is shown in the battery compartment.
- Remove spent batteries immediately from the device. There is an increased danger of leaks!
- Batteries must not be disposed of in the normal domestic waste!
- Every consumer is legally obliged to dispose of batteries in the proper manner!
- Keep batteries away from children; do not throw them in the fire, short-circuit them or take them apart.
- If the above instructions are not complied with, the batteries may discharge themselves beyond their end voltage. There is then a danger of

leaking. If the batteries have leaked inside your device, you should remove them immediately in order to prevent damage to the device!

- Always replace all batteries at the same time and only use the same type of batteries.

● Before getting started

Note: Press the count-up button **[8]** or the count-down-button **[9]** down to speed up the process of setting the values.

● Inserting / replacing the batteries

Note: When inserting the batteries, make sure that the contacts do not become bent.

- Open the battery compartment **[6]**, **[16]**, by pushing the battery compartment lid in the direction of the arrows.
- Take the battery compartment lid off.
- Insert the batteries in accordance with polarity.
- Close the battery compartment **[6]**, **[16]** by placing the battery compartment lid back on and pushing it in the opposite direction to the arrows until it perceptibly and audibly clicks into place.
- If the scales do not indicate any function, check that the batteries have been correctly inserted. You may have to reinsert them.

● Mounting the control unit

Note: You can mount the control unit onto the wall.

⚠ RISK OF ELECTRIC SHOCK! RISK OF DAMAGE TO PROPERTY! RISK OF INJURY! Make sure that you do not damage power cables, gas or water pipes when you drill into the wall. You may have to check with a cable locator before you drill into a wall.

- Place the scales about 10 cm away from the wall on a level and (scratch-) resistant surface.
- Select a place on the wall where you wish to mount the control unit.

Note: The control unit must not be mounted more than 1.5 m away from the scales. Otherwise no infrared connection can be established.

3. Mark the selected position with a pencil.
4. Use a drill to drill a hole in the wall. Please comply with the operating and safety instructions of the drill.
5. Place the dowel into the hole.
6. Screw the screw into the dowel. Make sure that the head of the screw sticks out of the wall by 5-7 mm.
7. Hang the control unit onto the wall using the mounting device [12].

● Preparing for use

● Setting the time

1. Press the time / confirmation button [7] down for about 3 seconds to be able to set the minutes.
Note: The hour display will flash.
2. Press the count-up button [8] or the count-down button [9] to set the hours.
3. Press the time / confirmation button [7] to confirm the entry and to be able to set the minutes.
Note: The minute display will flash.
4. Press the count-up button [8] or the count-down button [9] to set the minutes.
5. Press the time / confirmation button [7] to confirm your entry.

● Setting the alarm

1. Press the alarm-setting button [8] twice to activate the alarm and to display the alarm time.
Note: When the alarm is activated, the symbol  [18] will appear in the LC display [11].
2. Press the confirmation button [7] down for about 3 seconds to be able to set the hours.
Note: The hour display will flash.
3. Press the count-up button [8] or the count-down button [9] to set the hours.
4. Press the confirmation button [7] to confirm the entry and to be able to set the minutes.
Note: The minute display will flash.

5. Press the count-up button [8] or the count-down button [9] to set the minutes.
6. Press the confirmation button [7] to confirm the entry.
7. Wait for about 30 seconds or press the confirmation button [7] again in order to get back to time display.
8. Press the alarm-setting button [8] twice to deactivate the alarm.
Note: When the alarm is deactivated, the symbol  [18] in the LC display [11] will go out. The alarm time will go out.
9. Press the time / confirmation button [7] in order to get back to time display.

● Setting the unit of weight

Note: the legal weight unit is kg.

Note: You can choose to display the units of weight in either kilogrammes or "stones and pounds" or pound.

Note: When the unit of weight is changed, the unit of measurement for height also changes.

1. Pick the control unit up.
2. Make sure that the time is being displayed and press the button to set the units of weight / count-down button [9].
3. Press the time / confirmation button [7] down for about 3 seconds to be able to set the units of weight you wish to have.
Note: The pictograms of the units of weight [20], [21] will flash on the right-hand side of the LC display [11].
4. Press the count-up button [8] or the count-down button [9] to set the unit of weight.
5. Press the time / confirmation button [7] to confirm your entry.
6. Wait for about 30 seconds or press the time / confirmation button [7] again in order to get back to time display.

● Entering personal details

1. Pick the control unit up and make sure that the time is being displayed.

2. Press the time / confirmation button [7] in order to be able to enter your personal details.
Note: The letters "P-X" will flash in the LC display [11].
Note: "X" stands for the data profile of a person from the group of users.
3. Press the count-up button [8] or the count-down button [9], to select the data profile (1 - 12 or guest) you wish to have.
4. Press the time / confirmation button [7] to confirm the data profile you have selected and to set the unit of weight [20], [21].
5. Press the count-up button [8] or the count-down button [9] to set the unit of weight [20], [21].
6. Press the time / confirmation button [7] to confirm the entry and to be able to select your gender [34], [36].
7. Press the count-up button [8] or the count-down button [9] to select either male [36] or female [34].
8. Press the time / confirmation button [7] to confirm the entry and to be able to set your age.
9. Press the count-up button [8] or the count-down button [9] to set your age.
10. Press the time / confirmation button [7] to confirm the entry and to set your activity level [32], [33], [35], [37], [38].
11. Press the count-up button [8] or the count-down button [9] to set the appropriate activity level [32], [33], [35], [37], [38].
12. Press the time / confirmation button [7] to confirm the entry and to be able to set your height.
13. Press the count-up button [8] or the count-down button [9] to set your height.
14. Press the time / confirmation button [7] to confirm the entry.
15. Wait about 30 seconds to get back to time display.

● How to use

● Measuring weight

Note: Make sure that the control unit is at most 1.5 m from the scales when weighing yourself. Otherwise no infrared connection can be established. In this event, the LC display [11] will display "----".

- Place the scales down on a solid, flat horizontal surface.
 - Do not place the scales on a carpet. Otherwise this could influence the result. If necessary you can place a board on the carpet and place the scales upon this. A solid floor covering is required in order to achieve correct weight measurements.
1. Pick the control unit up and make sure that the time is being displayed.
 2. Press the count-down button [9] to be able to measure your weight only.
 3. Wait until "0,0 kg", "0,0 lb" or "0:0o st:lb" is displayed in the LC display [11].
 4. Step onto the platform [3] and stand still. Make sure that you are standing on both electrodes [2]. The weight measurement process will now start. The display will blink 3 times. The weight measured will be displayed and stored.
 5. Step down off the platform [3].
 6. Wait for about 30 seconds or press the time / confirmation button [7] in order to get back to time display.

● Measuring weight, calorie requirement, body fat, body water, bone and muscle percentages

Note: Make sure that the control unit is at most 1.5 m from the scales when weighing yourself. Otherwise no infrared connection can be established. In this event, the LC display [11] will display "----".
Tip: Always try to weigh yourself at the same time of day (best done in the morning), after going to the toilet, with an empty stomach and without any clothes on in order to get comparable results.

- Place the scales down on a solid, flat horizontal surface.
- Do not place the scales on a carpet. Otherwise this could influence the result. If necessary you can place a board on the carpet and place the scales upon this. A solid floor covering is required in order to achieve correct weight measurements.

1. Pick the control unit up and make sure that the time is being displayed.

Note: You can leave the control unit in the scales. Make sure that the metal contact  of the control unit lies precisely on the metal contact  of the scales.

2. Press the time / confirmation button  to be able to get the stored personal data.
3. Press the count-up button  or the count-down button  to get your profile memory.

Note: When your data profile is displayed, the LC display  will display all your personal settings.

4. Press the start button .
5. Wait until "0,0 kg", "0,0 lb" or "0:0o st:lb" is displayed in the LC display .
6. Step onto the platform  and stand still. Make sure that you are standing on both electrodes . The weight measurement process will now start. The display will blink 3 times. The weight measured will be displayed and stored. The LC display  will then display "00.0 000". After that, it will display the measured value three times in rotation in two groups.

Group 1: Calorie requirement  "kcal", body fat percentage   in % and body water percentage   in %.

Group 2: Weight "kg / lb / st:lb"  , muscle mass   in % and bone mass   in %.

After the values have been displayed, the time will automatically be displayed again.

● Retrieving last recorded values

Note: This function is only available for data profiles 1 - 12.

1. Make sure that the time is being displayed.
2. Press the start button  to be able to retrieve the stored readings.
3. Press the count-up button  or the count-down button  to select your data profile.
4. Press the time / confirmation button  and wait about 3-4 seconds to allow the last stored readings of the selected profile memory to be displayed. The readings will be displayed in rotation in two groups in the LC display  (see "Measuring weight, calorie requirement, body fat, body water, bone and muscle percentages", Point 6.).

● What the interpreted readings mean

LCD display					
Meaning	too thin	thin	normal	overweight	extremely overweight

● Trouble-shooting

Error display	Meaning	Solution
h-Lo	The control unit batteries are flat.	Change the batteries in line with the section called "Inserting / replacing the batteries".
b-Lo	The batteries in the scales are flat.	Change the batteries in line with the section called "Inserting / replacing the batteries".

Error display	Meaning	Solution
Err	The maximum load of 180 kg (400 lb) has been exceeded.	Step off the platform [3] of the scales in order to avoid injuries or damage to property.
ErrL	The body fat percentage measured is too low or is below the minimum measurable value of 4.0%.	-
Errh	The body fat percentage measured is too high or is above the maximum measurable value of 60.0%.	-
Various error messages appear in the LC display [11] .	Radio transmission devices (mobile phones, microwave ovens, etc.) interfere with the operation of the scale.	- Remove the radio transmission devices from proximity of the scale. - Remove the batteries from the battery compartment [6] and wait for about 10 seconds, then insert the battery back to the battery compartment [6] .

● **Cleaning**

- Clean the device using a slightly moist cloth. For more stubborn dirt you can also apply a weak soap solution to the cloth.
- Do not use chemical cleaners or scouring agents for cleaning the device.
- Never immerse the device in water or other liquids. Ignoring this advice may damage the device.

● **Disposal**

The packaging is made entirely of recyclable materials, which you may dispose of at local recycling facilities.

Contact your local refuse disposal authority for more details of how to dispose of your worn-out product.

To help protect the environment, please dispose of the product properly when it has reached the end of its useful life and not in the household waste. Information on collection points and their opening hours can be obtained from your local authority.

Faulty or used batteries must be recycled in accordance with Directive 2006 / 66 / EC. Please return the batteries and / or the device to the available collection points.

Environmental damage through incorrect disposal of the batteries!

Batteries may not be disposed of with the usual domestic waste. They may contain toxic heavy metals and are subject to hazardous waste treatment rules and regulations. The chemical symbols for heavy metals are as follows: Cd = cadmium, Hg = mercury, Pb = lead. That is why you should dispose of used rechargeable batteries at a local collection point.

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© by **ORFGEN Marketing**
Milomex Ltd.